• Welcome

Congratulations on your new iHealth Wireless Scale Lite with *Bluetooth*® 4.0 low-energy technology. This Scale tracks your weight and BMI and, along with the free companion “iHealth MyVitals” app, can help you stay motivated and on track for your weight goals.

• Mobile Device Compatibility

The iHealth Wireless Scale Lite works with the following devices with iOS version 5.0 or higher :

- iPhone 4S +
- iPad mini
- iPad (3rd generation) +
- iPod touch (5th generation) +

• Companion Mobile App

Scan the QR code and download the free “iHealth MyVitals” companion app or directly download the “iHealth MyVitals” app from the App Store. Follow the on-screen instructions to register and set up your iHealth ID.

iHealth MyVitals

• Access iHealth Cloud

Upon setting up your iHealth ID, you will also have access to a free, secure iHealth cloud account. Go to www.ihealthlabs.com, then click on “Sign In” to access your cloud account from PC or Mac using your iHealth ID.


• Install the Batteries

Open the battery compartment on the bottom of the Scale and install the batteries. Taking measurements twice a day, the battery power should last approximately three months.

• First Time Use

For your first weigh-in, follow these steps to ensure the Scale and the app are connected correctly and to synchronize the date and time stamp for your future weigh-ins.


1. Enable the *Bluetooth* on your mobile device (Settings->*Bluetooth*->On).

2. Tap the Scale to turn it on. The *Bluetooth* icon on the Scale will begin to flash.

3. Launch the app to initiate the connection. The *Bluetooth* icon will light up and stop flashing when a successful connection is established. The Scale is now connected and ready for your weigh-in.

• Your Next Weigh-in

Subsequent daily weigh-ins can be taken without being connected to your mobile device, as your weight data is stored in the Scale and uploaded to the app upon your next connection. Real-time weight results are viewable in the app if the Scale is connected to your mobile device and the app is running when your weigh-in occurs.

• Software updates

From time to time, the Scale’s software is updated to enhance performance. When a software update is available, a notification message will appear in the app. Complete the following steps to install the update:

- ❶ Follow steps 1-3 in the “First-time Use” section.
- ❷ You will see a notification message in the app asking for permission to proceed. Select “Yes”.
- ❸ The update will begin to download. You will see a cursor on the scale display moving from left to right. The scale will display “0.0” when the update is complete.

Visit www.ihealthlabs.com to obtain additional product information. For Customer Service, please call +1(855)816-7705.

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