

iHealth™

Wireless Scale (HS3)

QUICK START GUIDE

DOWNLOAD THE FREE iHealth MyVitals

Prior to first use, download and install the “iHealth MyVitals” from the App Store.

PRE-SET UP INSTRUCTIONS


- iOS device is compatible and is version V5.0 or higher.

INSTALL BATTERIES

Open battery cover on the back of iHealth Wireless Scale and install four AAA batteries.




SET UP THE *BLUETOOTH* CONNECTION

STEP 1: Turn on the Scale by stepping on it until “0.0” is displayed.



STEP 2: Turn the *Bluetooth* setting “On” on your mobile device.

STEP 3: Wait until the model name, “iHealth HS3XXXXX”, and “Not Paired” appear on the device screen. Select the model name “iHealth HS3XXXXX” to pair and connect. (Note: it may take up to 30 seconds for the device to detect the *Bluetooth* signal).



STEP 4: After the *Bluetooth* connection is successful, launch the App.



MEASUREMENT PROCEDURES

Step on iHealth Wireless Scale. Once the Scale detects that there is no more movement, weight information will display on the App screen.



Visit www.ihealthlabs.com to obtain additional product information. For Customer Services, please call +1(855)816-7705.